

OVERVIEW:

The Sys*Stim 206 is a one-channel neuromuscular stimulator that features four different waveforms in a light-weight package. The narrow pulse is suited for TENS while the wide and AC pulses can be used for muscle stimulation. Direct Current can be used to stimulate denervated muscles or perform iontophoresis. The Sys*Stim 206 is an ideal companion to the Sonicator 716 for combination therapy.

The Mettler Sys*Stim 206 has four preset stimulation modes plus Surge and Reverse Polarity.

Each mode allows the most beneficial waveform for different treatments
Each mode is balanced for maximum patient comfort
The Sys*Stim 206 overcomes current breakdown caused by skin resistance
The waveform you select is always the one you get
Portable, so you can take it anywhere
LED display indicate output intensity in peak volts or milliamps
LED displays treatment time remaining