

PRODUCT DESCRIPTION:

The Sonicator 716 is a lightweight, portable ultrasound that features a 10 cm² / 1 MHz applicator. It is ideal for treating shoulders, backs and other large, muscular areas. The Sonicator 716 is built to meet all international safety standards, so it can travel anywhere and deliver world-class quality and performance.

Weight: 5.1 pounds (2.3 kg.)

Dimensions: 4.3 in (H) x 6 in (D) x 13.4 in (L), (11 cm (H) x 15 cm (D) x 34 cm (L)

Mettler Ultrasound -Direct Crystal-to-Patient

All Mettler Ultrasound - Sonicator ultrasound applicators feature unique crystal-to-patient technology that delivers ultrasound directly to the patient without a metal interface as well as watertight characteristics for underwater therapy.